

Fitzpatrick Skin Type	Description	Starting Fluences (J) (Test Spots)	Max J	Pulse duration
I	Light white skin always burns, never tans	26J/30 28J/30 30J/30 32J/30 34J/30	34J	Auto (finer hair) or 30 ms
II	Light white skin always burns, sometimes tans	26J/30 28J/30 30J/30 32J/30 34J/30	34J	Auto (finer hair) or 30ms-high
III	Medium white skin sometimes burns, always tans	24J/30 26J/30 28J/30 30J/30	30J	30ms
IV	Dark/olive white and Asian skin rarely burns, always tans	Extremely Careful 24 J/100 25J/100 26J/100 27J/100 28/100 max	28J	100-400ms
V	Light brown skin	Extremely Careful 20 J/400 or 10J/100 No Higher		400 ms
VI	Medium to dark brown, African and African-American skin	We don't treat		