

Laser Facial - Yag Laser - Lyra I - this is to defuse redness (browns is IPL, combo brown and red is IPL)

Skin Type 1-3	5mm (3-5 passes)	13-17J - 30 pulse
Skin Type 4	5 mm (2-3 passes)	13J - 100 pulse

Spider Veins - Yag Laser - Lyra I - If blanching (white spots) reduce power, always double cool - gel

Skin type 1-3 only		250-370J increase 10% in needed - single pass
Legs	3mm (1-2 mm vessel)	150-200J
Legs	5mm (2-4mm vessel)	180J-200J - single pass only
face	1.5mm	100-180J (no seam of nose) single pass only, avoid pitting.
nose	1.5 mm	Follow with 19J 5mm 2pps 5 passes

also refer to Lyra settings sheet for 1.5 mm spot
always be careful especially on nose, reduce power on nose, do not treat crease of nose

Yag - Hair Removal Lyra-I

Skin Type 4	30-50J	100 pulse	use gel and air go slow
Skin type 5	26-35 J	100 pulse	use gel and air go slow
Skin type 6	We do not treat		